

ABX Minocal



LOT

CX 432

CAL



(Exp.) 2019-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.6	9.0	9.1	9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.35	4.28	4.37	4.34	4.38					± 0.06
HB HGB	g/dl	13.3	13.3	13.2	13.2	13.3					± 0.2
	g/l	133	133	132	132	133					± 2
	mmol/l	8.26	8.26	8.20	8.20	8.26					± 0.12
HT HCT	%	37.3	37.2	37.6	37.8	37.9					± 1.0
	l/l	0.373	0.372	0.376	0.378	0.379					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	245	239	249	244					± 10
VPM MPV	µm ³ ; fl	8.8	8.6	N/A	9.2	8.0					± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal



LOT

CX 432
Rev 1

CAL

(Exp.) **2019-04-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1		9.2	9.2	9.4	9.0	9.4	9.1	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.35		4.32	4.36	4.39	4.37	4.39	4.36	4.36	± 0.06
HB HGB	g/dl	13.1		13.2	13.2	13.3	13.2	13.3	13.2	13.2	± 0.2
	g/l	131		132	132	133	132	133	132	132	± 2
	mmol/l	8.14		8.20	8.20	8.26	8.20	8.26	8.20	8.20	± 0.12
HT HCT	%	37.2		35.4	36.2	36.7	35.8	36.4	36.2	36.2	± 1.0
	l/l	0.372		0.354	0.362	0.367	0.358	0.364	0.362	0.362	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245		243	240	250	252	250	251	249	± 10
VPM MPV	µm ³ ; fl	9.0		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

1300060880
BACK / VERSO
Ref: TEMP-0387 Rev.46

ABX Minocal



LOT

CX 432

CAL



(Exp.) 2019-04-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40								± 0.06
HB HGB	g/dl	13.3								± 0.2
	g/l	133								± 2
	mmol/l	8.26								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	253								± 10
VPM MPV	µm ³ ; fl	N/A								N/A
VGM MCV	µm ³ ; fl	83.5								± 4.0

1300060884
Ref: TEMP-0387 Rev.46