

ABX Minocal



(Exp.) 2020-11-05

(YYYY - MM - DD)

LOT		CX 451		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.2	8.9	9.2	9.3						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.47	4.46	4.52	4.51						± 0.06
HB HGB	g/dl	13.8	14.0	13.8	14.0	13.9						± 0.2
	g/l	138	140	138	140	139						± 2
	mmol/l	8.57	8.69	8.57	8.69	8.63						± 0.12
HT HCT	%	38.8	39.6	39.2	39.6	39.7						± 1.0
	l/l	0.388	0.396	0.392	0.396	0.397						± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	242	238	246	245						± 10
VPM MPV	µm ³ ; fl	7.8	8.0	N/A	8.6	7.4						± 0.5

Ref. TEMP-0387 Rev.46 FRONT / RECTO 1300060880

LOT		CX 451		CAL		(Exp.) 2020-11-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4		9.4	9.1	9.7	9.5	9.6	9.3	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51		4.43	4.50	4.51	4.50	4.51	4.50	4.50	± 0.06
HB HGB	g/dl	13.7		13.8	13.9	13.9	13.8	13.9	13.9	13.9	± 0.2
	g/l	137		138	139	139	138	139	139	139	± 2
	mmol/l	8.51		8.57	8.63	8.63	8.57	8.63	8.63	8.63	± 0.12
HT HCT	%	39.2		37.0	38.0	37.9	37.4	37.9	38.0	38.3	± 1.0
	l/l	0.392		0.370	0.380	0.379	0.374	0.379	0.380	0.383	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		238	240	243	243	243	247	243	± 10
VPM MPV	µm ³ ; fl	8.3		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300060880
 BACK / VERSO
 Ref. TEMP-0387 Rev.46

ABX Minocal

LOT		CX 451		CAL		(Exp.) 2020-11-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
			YUMIZEN								
			H1500 H2500								
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90								± 0.20
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50								± 0.06
HB	HGB	g/dl	13.8								± 0.2
		g/l	138								± 2
		mmol/l	8.57								± 0.12
HT	HCT	%	N/A								N/A
		l/l	N/A								N/A
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	240								± 10
VPM	MPV	µm ³ ; fl	N/A								N/A
VGM	MCV	µm ³ ; fl	85.0								± 4.0

1300060884

Ref: TEMP-0387 Rev.46

 1300060884
 Ref: TEMP-0387 Rev.46