

ABX Minocal

LOT		CX 437		CAL		(Exp.) 2019-09-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.5	8.8	8.8	8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.48	4.45	4.45	4.46					± 0.06
HB HGB	g/dl	13.5	13.7	13.5	13.6	13.6					± 0.2
	g/l	135	137	135	136	136					± 2
	mmol/l	8.38	8.51	8.38	8.45	8.45					± 0.12
HT HCT	%	38.3	38.8	38.5	38.7	38.6					± 1.0
	l/l	0.383	0.388	0.385	0.387	0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	242	235	247	244					± 10
VPM MPV	µm ³ ; fl	8.0	8.0	N/A	8.6	7.6					± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal


 (Exp.) 2019-09-05
 (YYYY - MM - DD)

LOT		CX 437		CAL		(Exp.) 2019-09-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9		9.1	8.9	9.3	8.9	9.3	8.9	8.9	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45		4.40	4.46	4.49	4.44	4.49	4.47	4.46	± 0.06	
HB HGB	g/dl	13.4		13.5	13.7	13.5	13.5	13.5	13.6	13.7	± 0.2	
	g/l	134		135	137	135	135	135	136	137	± 2	
	mmol/l	8.32		8.38	8.51	8.38	8.38	8.38	8.45	8.51	± 0.12	
HT HCT	%	38.3		36.1	37.5	37.3	36.6	37.3	37.3	37.5	± 1.0	
	l/l	0.383		0.361	0.375	0.373	0.366	0.373	0.373	0.375	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242		233	240	248	243	248	240	240	± 10	
VPM MPV	µm ³ ; fl	8.5		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300060880
 BACK / VERSO
 Ref: TEMP-0387 Rev.46

ABX Minocal


LOT
CX 437
CAL

(Exp.) 2019-09-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51									± 0.06
HB HGB	g/dl	13.5									± 0.2
	g/l	135									± 2
	mmol/l	8.38									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	83.8									± 4.0

1300060884

Ref: TEMP-0387 Rev.46