

ABX Minocal



(Exp.) 2019-08-05

(YYYY - MM - DD)

LOT		CX 436		CAL		(Exp.) 2019-08-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.4	8.7	8.7	8.8					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.55	4.54	4.54	4.57					± 0.06		
HB HGB	g/dl	13.5	13.8	13.6	13.7	13.8					± 0.2		
	g/l	135	138	136	137	138					± 2		
	mmol/l	8.38	8.57	8.45	8.51	8.57					± 0.12		
HT HCT	%	38.4	38.7	38.6	38.8	38.8					± 1.0		
	l/l	0.384	0.387	0.386	0.388	0.388					± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	237	230	238	236					± 10		
VPM MPV	µm ³ ; fl	8.7	8.4	N/A	9.2	8.0					± 0.5		

 1300060880
 FRONT / RECTO
 Ref. TEMP-0387 Rev.46

ABX Minocal


 (Exp.) 2019-08-05
 (YYYY - MM - DD)

LOT		CX 436		CAL		(Exp.) 2019-08-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8		8.8	8.8	9.2	8.8	9.2	8.8	8.8	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53		4.53	4.58	4.60	4.55	4.60	4.57	4.58	± 0.06	
HB HGB	g/dl	13.5		13.6	13.8	13.7	13.6	13.7	13.7	13.8	± 0.2	
	g/l	135		136	138	137	136	137	137	138	± 2	
	mmol/l	8.38		8.45	8.57	8.51	8.45	8.51	8.51	8.57	± 0.12	
HT HCT	%	38.1		36.9	37.6	37.7	36.9	37.5	37.3	37.6	± 1.0	
	l/l	0.381		0.369	0.376	0.377	0.369	0.375	0.373	0.376	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		227	226	240	235	240	235	233	± 10	
VPM MPV	µm ³ ; fl	9.0		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300060880
 BACK / VERSO
 Ref. TEMP-0387 Rev.46

ABX Minocal



LOT		CX 436		CAL		(Exp.) 2019-08-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.55									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63									± 0.06	
HB	HGB	g/dl	13.7									± 0.2	
		g/l	137									± 2	
		mmol/l	8.51									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	246									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	82.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.46