

ABX Minocal



(Exp.) 2017-12-05

(YYYY - MM - DD)

LOT		CX 416		CAL		(Exp.) 2017-12-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2	9.1	9.1	9.2					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.80	4.79	4.82	4.78	4.79					± 0.06		
HB HGB	g/dl	13.8	13.9	13.9	13.8	13.9					± 0.2		
	g/l	138	139	139	138	139					± 2		
	mmol/l	8.57	8.63	8.63	8.57	8.63					± 0.12		
HT HCT	%	39.8	40.0	40.0	39.7	39.8					± 1.0		
	l/l	0.398	0.400	0.400	0.397	0.398					± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	251	252	259	257					± 10		
VPM MPV	µm ³ ; fl	9.3	9.2	N/A	9.3	8.7					± 0.5		

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1219930090 -A

ABX Minocal


 (Exp.) 2017-12-05
 (YYYY - MM - DD)

LOT		CX 416		CAL		(Exp.) 2017-12-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3		9.4	9.2	9.6	9.3	9.6	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.80		4.77	4.80	4.84	4.81	4.84	4.80	4.80	± 0.06
HB HGB	g/dl	13.8		13.7	13.8	13.8	14.0	13.8	13.8	13.8	± 0.2
	g/l	138		137	138	138	140	138	138	138	± 2
	mmol/l	8.57		8.51	8.57	8.57	8.69	8.57	8.57	8.57	± 0.12
HT HCT	%	39.8		37.3	37.9	39.0	38.0	38.7	38.1	38.2	± 1.0
	l/l	0.398		0.373	0.379	0.390	0.380	0.387	0.381	0.382	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	265		245	250	258	262	258	260	260	± 10
VPM MPV	µm ³ ; fl	9.4		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1219930090 -A
 BACK / VERSO
 Ref: TEMP-0387 Rev.44

ABX Minocal

LOT		CX 416		CAL		(Exp.) 2017-12-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
			YUMIZEN								
			H1500 H2500								
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9								± 0.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.87								± 0.06
HB	HGB	g/dl	13.8								± 0.2
		g/l	138								± 2
		mmol/l	8.57								± 0.12
HT	HCT	%	N/A								N/A
		l/l	N/A								N/A
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	265								± 10
VPM	MPV	µm ³ ; fl	N/A								N/A
VGM	MCV	µm ³ ; fl	83								± 4

 1300032517-A
 Ref: TEMP-0387 Rev.44