

ABX Minocal



(Exp.) 2019-07-05

(YYYY - MM - DD)

LOT		CX 435		CAL		(Exp.) 2019-07-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.4	8.7	8.8	9.0					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.55	4.54	4.55	4.58					± 0.06		
HB HGB	g/dl	13.6	13.5	13.4	13.5	13.6					± 0.2		
	g/l	136	135	134	135	136					± 2		
	mmol/l	8.45	8.38	8.32	8.38	8.45					± 0.12		
HT HCT	%	38.3	38.2	38.1	38.4	38.5					± 1.0		
	l/l	0.383	0.382	0.381	0.384	0.385					± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	242	234	242	241					± 10		
VPM MPV	µm ³ ; fl	8.2	7.9	N/A	8.6	7.5					± 0.5		

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal


 (Exp.) 2019-07-05
 (YYYY - MM - DD)

LOT		CX 435		CAL		(Exp.) 2019-07-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9		8.9	8.9	9.2	8.8	9.2	8.9	8.9	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54		4.51	4.58	4.60	4.57	4.60	4.58	4.58	± 0.06
HB HGB	g/dl	13.3		13.4	13.5	13.5	13.4	13.5	13.5	13.6	± 0.2
	g/l	133		134	135	135	134	135	135	136	± 2
	mmol/l	8.26		8.32	8.38	8.38	8.32	8.38	8.38	8.45	± 0.12
HT HCT	%	37.8		36.1	36.9	37.4	36.3	37.0	36.8	37.0	± 1.0
	l/l	0.378		0.361	0.369	0.374	0.363	0.370	0.368	0.370	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		240	236	243	247	243	244	240	± 10
VPM MPV	µm ³ ; fl	8.2		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300060880
 BACK / VERSO
 Ref: TEMP-0387 Rev.46

ABX Minocal


LOT
CX 435
CAL

(Exp.) 2019-07-05

 Rev **1**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62									± 0.06
HB HGB	g/dl	13.4									± 0.2
	g/l	134									± 2
	mmol/l	8.32									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	81.0									± 4.0

1300060884

Ref: TEMP-0387 Rev.46