

ABX Minocal

LOT
CX 464
CAL

(Exp.) 2021-12-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2			9.0				± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.55			4.50				± 0.06
HB HGB	g/dl	13.5	13.6			13.4				± 0.2
	g/l	135	136			134				± 2
	mmol/l	8.38	8.45			8.32				± 0.12
HT HCT	%	38.6	38.7			38.5				± 1.0
	l/l	0.386	0.387			0.385				± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	248			245				± 10
VPM MPV	µm ³ ; fl	8.8	8.6			8.1				± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300090308

ABX Minocal


 (Exp.) 2021-12-05
 (YYYY - MM - DD)

LOT		CX 464		CAL		(Exp.) 2021-12-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0	9.1	9.0	9.4	9.0	9.2	9.0	9.0	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.45	4.45	4.48	4.48	4.48	4.51	4.52	4.48	± 0.06	
HB HGB	g/dl	13.3	13.2	13.3	13.4	13.3	13.3	13.4	13.4	13.4	± 0.2	
	g/l	133	132	133	134	133	133	134	134	134	± 2	
	mmol/l	8.26	8.20	8.26	8.32	8.26	8.26	8.32	8.32	8.32	± 0.12	
HT HCT	%	37.7	37.4	36.3	36.5	36.3	36.1	36.5	36.6	36.5	± 1.0	
	l/l	0.377	0.374	0.363	0.365	0.363	0.361	0.365	0.366	0.365	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	249	248	244	240	246	246	250	245	240	± 10	
VPM MPV	µm ³ ; fl	9.1	9.0	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref. TEMP-0387 Rev.47

ABX Minocal



LOT		CX 464		CAL		(Exp.) 2021-12-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52									± 0.06	
HB	HGB	g/dl	13.4									± 0.2	
		g/l	134									± 2	
		mmol/l	8.32									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	254									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	82.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.47