

# ABX Minocal



**LOT** CX 102  
Rev 1

**CAL**

(Exp.) 2012-11-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	$10^3/mm^3 ; 10^9/l$	9.8	9.4	9.5	9.6	9.6		10.0	9.5	9.4	± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.70	4.70	4.70	4.72	4.72		4.74	4.70	4.70	± 0.06
HB HGB	g/dl	13.8	13.8	13.7	13.7	13.8		13.7	13.6	13.7	± 0.2
	g/l	138	138	137	137	138		137	136	137	± 2
	mmol/l	8.57	8.57	8.51	8.51	8.57		8.51	8.45	8.51	± 0.12
HT HCT	%	39.2	39.3	39.5	39.6	39.6		38.2	37.6	38.1	± 1.0
	l/l	0.392	0.393	0.395	0.396	0.396		0.382	0.376	0.381	± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	257	247	259	255	255		252	258	258	± 10
VPM MPV	$\mu m^3 ; fl$	7.8	7.8	N/A	7.8	7.7		N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.22 FRONT / RECTO 9930062-A

# ABX Minocal



**LOT** CX 102  
Rev 1

**CAL**

(Exp.) 2012-11-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMICRP	CRP CRP 200	8P	18P	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	$10^3/mm^3 ; 10^9/l$	<b>9.4</b>				<b>10.0</b>	<b>9.4</b>	<b>10.0</b>	<b>9.5</b>	<b>9.4</b>	± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	<b>4.69</b>				<b>4.74</b>	<b>4.71</b>	<b>4.74</b>	<b>4.70</b>	<b>4.70</b>	± 0.06
HB HGB	g/dl	<b>13.6</b>				<b>13.7</b>	<b>13.7</b>	<b>13.7</b>	<b>13.6</b>	<b>13.7</b>	± 0.2
	g/l	<b>136</b>				<b>137</b>	<b>137</b>	<b>137</b>	<b>136</b>	<b>137</b>	± 2
	mmol/l	<b>8.45</b>				<b>8.51</b>	<b>8.51</b>	<b>8.51</b>	<b>8.45</b>	<b>8.51</b>	± 0.12
HT HCT	%	<b>39.6</b>				<b>38.2</b>	<b>36.7</b>	<b>38.2</b>	<b>37.6</b>	<b>38.1</b>	± 1.0
	l/l	<b>0.396</b>				<b>0.382</b>	<b>0.367</b>	<b>0.382</b>	<b>0.376</b>	<b>0.381</b>	± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	<b>245</b>				<b>252</b>	<b>254</b>	<b>252</b>	<b>258</b>	<b>258</b>	± 10
VPM MPV	$\mu m^3 ; fl$	<b>8.2</b>				N/A	N/A	N/A	N/A	N/A	± 0.5

9930062-A  
BACK / VERSO  
Ref: TEMP-0387 Rev.22