

ABX Minocal


LOT
CX 033
CAL

(Exp.) 2013-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	9.3	9.0	9.1	9.0	9.2		9.4	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.57	4.59	4.56	4.56	4.56		4.60	4.57	4.56	± 0.06
HB HGB	g/dl	13.7	13.6	13.5	13.6	13.5		13.6	13.5	13.6	± 0.2
	g/l	137	136	135	136	135		136	135	136	± 2
	mmol/l	8.51	8.45	8.38	8.45	8.38		8.45	8.38	8.45	± 0.12
HT HCT	%	39.0	39.2	39.0	39.0	38.9		37.7	37.4	37.6	± 1.0
	l/l	0.390	0.392	0.390	0.390	0.389		0.377	0.374	0.376	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	240	233	238	244	236		240	238	236	± 10
VPM MPV	µm ³ ; fl	8.3	8.3	N/A	8.6	8.0		N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.31 FRONT / RECTO 9930077 -A

ABX Minocal


LOT
CX 033
CAL

(Exp.) 2013-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	CRP CRP 200	8P	18P	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	8.7				9.4	9.1	9.4	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.58				4.60	4.55	4.60	4.57	4.56	± 0.06
HB HGB	g/dl	13.3				13.6	13.6	13.6	13.5	13.6	± 0.2
	g/l	133				136	136	136	135	136	± 2
	mmol/l	8.26				8.45	8.45	8.45	8.38	8.45	± 0.12
HT HCT	%	39.3				37.7	36.2	37.7	37.4	37.6	± 1.0
	l/l	0.393				0.377	0.362	0.377	0.374	0.376	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	232				240	238	240	238	236	± 10
VPM MPV	µm ³ ; fl	8.4				N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.31 BACK / VERSO 9930077 -A