

ABX Minocal



(Exp.) 2018-01-05

(YYYY - MM - DD)

LOT		CX 417		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.3	9.0	9.3	9.3						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.67	4.63	4.70	4.63	4.65						± 0.06
HB HGB	g/dl	13.7	13.7	13.7	13.6	13.7						± 0.2
	g/l	137	137	137	136	137						± 2
	mmol/l	8.51	8.51	8.51	8.45	8.51						± 0.12
HT HCT	%	39.0	38.7	39.2	38.7	39.0						± 1.0
	l/l	0.390	0.387	0.392	0.387	0.390						± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	254	251	253	253	253						± 10
VPM MPV	µm ³ ; fl	8.8	8.6	N/A	8.8	8.1						± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1219930090 -A

ABX Minocal


 (Exp.) 2018-01-05
 (YYYY - MM - DD)

LOT		CX 417		CAL		(Exp.) 2018-01-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3		9.6	9.3	9.5	9.3	9.5	9.3	9.3	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.65		4.66	4.66	4.70	4.72	4.70	4.67	4.68	± 0.06
HB HGB	g/dl	13.5		13.6	13.6	13.7	13.8	13.7	13.6	13.6	± 0.2
	g/l	135		136	136	137	138	137	136	136	± 2
	mmol/l	8.38		8.45	8.45	8.51	8.57	8.51	8.45	8.45	± 0.12
HT HCT	%	38.5		36.8	37.0	38.1	37.3	37.8	37.1	37.4	± 1.0
	l/l	0.385		0.368	0.370	0.381	0.373	0.378	0.371	0.374	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	257		250	250	255	261	255	258	255	± 10
VPM MPV	µm ³ ; fl	8.7		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1219930090 -A

ABX Minocal

LOT
CX 417
CAL

(Exp.) 2018-01-05

 Rev **1**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.74									± 0.06
HB HGB	g/dl	13.6									± 0.2
	g/l	136									± 2
	mmol/l	8.45									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	270									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	79									± 4

 1300032517-A
 Ref: TEMP-0387 Rev.44