

ABX Minocal



(Exp.) 2022-07-05

(YYYY - MM - DD)

LOT		CX 471		CAL		(Exp.) 2022-07-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.8			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.49			4.52					± 0.06
HB HGB	g/dl	13.4	13.6			13.4					± 0.2
	g/l	134	136			134					± 2
	mmol/l	8.32	8.45			8.32					± 0.12
HT HCT	%	38.1	38.6			38.6					± 1.0
	l/l	0.381	0.386			0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	251	244			251					± 10
VPM MPV	µm ³ ; fl	8.8	8.6			7.8					± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.47

ABX Minocal


 (Exp.) 2022-07-05
 (YYYY - MM - DD)

LOT		CX 471				CAL		(Exp.) 2022-07-05 (YYYY - MM - DD)			
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.2	9.0	9.3	9.1	9.3	9.1	9.0	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.45	4.40	4.49	4.52	4.49	4.52	4.51	4.49	± 0.06
HB HGB	g/dl	13.2	13.3	13.2	13.5	13.4	13.3	13.4	13.4	13.5	± 0.2
	g/l	132	133	132	135	134	133	134	134	135	± 2
	mmol/l	8.20	8.26	8.20	8.38	8.32	8.26	8.32	8.32	8.38	± 0.12
HT HCT	%	37.6	37.6	36.5	36.8	36.8	36.4	36.8	37.0	36.8	± 1.0
	l/l	0.376	0.376	0.365	0.368	0.368	0.364	0.368	0.370	0.368	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	253	242	245	250	247	250	245	245	± 10
VPM MPV	µm ³ ; fl	8.9	8.9	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300090308

ABX Minocal



LOT		CX 471		CAL		(Exp.) 2022-07-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53									± 0.06	
HB	HGB	g/dl	13.5									± 0.2	
		g/l	135									± 2	
		mmol/l	8.38									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	252									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	84.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.47