

ABX Minocal



(Exp.) 2021-01-05

(YYYY - MM - DD)

LOT		CX 453		CAL		(Exp.) 2021-01-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0	8.9	9.0	9.1					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.55	4.55	4.53	4.50					± 0.06		
HB HGB	g/dl	13.9	14.2	14.1	14.0	14.0					± 0.2		
	g/l	139	142	141	140	140					± 2		
	mmol/l	8.63	8.82	8.76	8.69	8.69					± 0.12		
HT HCT	%	39.8	40.0	40.3	40.1	39.8					± 1.0		
	l/l	0.398	0.400	0.403	0.401	0.398					± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	240	235	240	240					± 10		
VPM MPV	µm ³ ; fl	8.2	8.2	N/A	9.5	7.8					± 0.5		

 1300060880
 FRONT / RECTO
 Ref. TEMP-0387 Rev.46

ABX Minocal


 (Exp.) 2021-01-05
 (YYYY - MM - DD)

LOT		CX 453		CAL		(Exp.) 2021-01-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3		9.3	9.1	9.4	9.2	9.3	9.2	9.1	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52		4.49	4.53	4.55	4.51	4.55	4.52	4.53	± 0.06
HB HGB	g/dl	13.8		13.9	14.1	14.0	13.9	14.0	14.1	14.1	± 0.2
	g/l	138		139	141	140	139	140	141	141	± 2
	mmol/l	8.57		8.63	8.76	8.69	8.63	8.69	8.76	8.76	± 0.12
HT HCT	%	39.1		37.9	38.5	38.4	37.7	38.4	38.4	38.7	± 1.0
	l/l	0.391		0.379	0.385	0.384	0.377	0.384	0.384	0.387	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237		234	238	244	239	244	243	242	± 10
VPM MPV	µm ³ ; fl	8.5		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300060880
 BACK / VERSO
 Ref. TEMP-0387 Rev.46

ABX Minocal

LOT		CX 453		CAL		(Exp.) 2021-01-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56									± 0.06	
HB	HGB	g/dl	14.0									± 0.2	
		g/l	140									± 2	
		mmol/l	8.69									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	242									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	85.5									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.46