

ABX Minocal

LOT		CX 413		CAL		(Exp.) 2017-09-05 (YYYY - MM - DD)					
		Rev 2									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1	9.0	8.9	9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.70	4.68	4.65	4.63	4.69					± 0.06
HB HGB	g/dl	13.9	13.8	13.8	13.7	13.9					± 0.2
	g/l	139	138	138	137	139					± 2
	mmol/l	8.63	8.57	8.57	8.51	8.63					± 0.12
HT HCT	%	39.7	39.8	39.5	39.4	39.9					± 1.0
	l/l	0.397	0.398	0.395	0.394	0.399					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	251	252	250	255	256					± 10
VPM MPV	µm ³ ; fl	8.8	8.4	N/A	8.9	8.1					± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1219930090 -A

ABX Minocal

LOT		CX 413		CAL		(Exp.) 2017-09-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2		9.3	9.0	9.3	9.2	9.3	9.1	9.0	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.66		4.63	4.70	4.71	4.67	4.71	4.67	4.70	± 0.06
HB HGB	g/dl	13.8		13.8	14.0	13.8	13.9	13.8	13.8	13.9	± 0.2
	g/l	138		138	140	138	139	138	138	139	± 2
	mmol/l	8.57		8.57	8.69	8.57	8.63	8.57	8.57	8.63	± 0.12
HT HCT	%	39.1		37.2	38.3	38.4	37.6	38.4	38.1	38.3	± 1.0
	l/l	0.391		0.372	0.383	0.384	0.376	0.384	0.381	0.383	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	263		246	255	257	258	257	254	255	± 10
VPM MPV	µm ³ ; fl	8.7		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1219930090 -A

ABX Minocal

LOT		CX 413		CAL		(Exp.) 2017-09-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.70									± 0.06
HB HGB	g/dl	13.7									± 0.2
	g/l	137									± 2
	mmol/l	8.51									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	260									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	83									± 4

 1300032517-A
 Ref: TEMP-0387 Rev.44