

ABX Minocal


LOT
CX 452
CAL

(Exp.) 2020-12-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.0	9.0	9.2				± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.55	4.50	4.49	4.54				± 0.06
HB HGB	g/dl	13.8	14.1	13.8	13.8	13.9				± 0.2
	g/l	138	141	138	138	139				± 2
	mmol/l	8.57	8.76	8.57	8.57	8.63				± 0.12
HT HCT	%	39.5	40.0	39.6	39.5	39.5				± 1.0
	l/l	0.395	0.400	0.396	0.395	0.395				± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	237	232	240	237				± 10
VPM MPV	µm ³ ; fl	8.0	7.9	N/A	8.6	7.4				± 0.5

Ref. TEMP-0387 Rev.46 FRONT / RECTO 1300060880

LOT
CX 452
CAL

(Exp.) 2020-12-05
 (YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4		9.4	9.1	9.6	9.3	9.5	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55		4.50	4.52	4.55	4.52	4.55	4.55	4.54	± 0.06
HB HGB	g/dl	13.8		13.8	13.9	13.9	13.8	13.9	13.9	13.9	± 0.2
	g/l	138		138	139	139	138	139	139	139	± 2
	mmol/l	8.57		8.57	8.63	8.63	8.57	8.63	8.63	8.63	± 0.12
HT HCT	%	39.1		37.8	38.2	38.2	37.7	38.2	38.4	38.4	± 1.0
	l/l	0.391		0.378	0.382	0.382	0.377	0.382	0.384	0.384	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237		240	236	240	242	240	244	245	± 10
VPM MPV	µm ³ ; fl	8.3		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060880

ABX Minocal



LOT		CX 452		CAL		(Exp.) 2020-12-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56									± 0.06	
HB	HGB	g/dl	13.9									± 0.2	
		g/l	139									± 2	
		mmol/l	8.63									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	240									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	85.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.46