

ABX Minocal



LOT		CX 451		CAL		(Exp.) 2020-11-05 (YYYY - MM - DD)					
Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.2	8.9	9.2	9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.47	4.46	4.52	4.51					± 0.06
HB HGB	g/dl	13.8	14.0	13.8	14.0	13.9					± 0.2
	g/l	138	140	138	140	139					± 2
	mmol/l	8.57	8.69	8.57	8.69	8.63					± 0.12
HT HCT	%	38.8	39.6	39.2	39.6	39.7					± 1.0
	l/l	0.388	0.396	0.392	0.396	0.397					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	242	238	246	245					± 10
VPM MPV	µm ³ ; fl	7.8	8.0	N/A	8.6	7.4					± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

 1300060880
 FRONT / RECTO
 Ref. TEMP-0387 Rev.46

LOT		CX 451		CAL		(Exp.) 2020-11-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4		9.4	9.1	9.7	9.5	9.6	9.3	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51		4.43	4.50	4.51	4.50	4.51	4.50	4.50	± 0.06
HB HGB	g/dl	13.7		13.8	13.9	13.9	13.8	13.9	13.9	13.9	± 0.2
	g/l	137		138	139	139	138	139	139	139	± 2
	mmol/l	8.51		8.57	8.63	8.63	8.57	8.63	8.63	8.63	± 0.12
HT HCT	%	39.2		37.0	38.0	37.9	37.4	37.9	38.0	38.3	± 1.0
	l/l	0.392		0.370	0.380	0.379	0.374	0.379	0.380	0.383	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		238	240	243	243	243	247	243	± 10
VPM MPV	µm ³ ; fl	8.3		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300060880
 BACK / VERSO
 Ref. TEMP-0387 Rev.46

ABX Minocal



LOT		CX 451		CAL		(Exp.) 2020-11-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50									± 0.06
HB HGB	g/dl	13.8									± 0.2
	g/l	138									± 2
	mmol/l	8.57									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	85.0									± 4.0

 1300060884
 Ref: TEMP-0387 Rev.46