

ABX Minocal


 (Exp.) 2019-03-05
 (YYYY - MM - DD)

LOT		CX 431		CAL		(Exp.) 2019-03-05 (YYYY - MM - DD)							
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.5	8.6	8.8	8.9					± 0.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.46	4.50	4.50	4.53					± 0.06		
HB HGB	g/dl	13.4	13.4	13.3	13.4	13.4					± 0.2		
	g/l	134	134	133	134	134					± 2		
	mmol/l	8.32	8.32	8.26	8.32	8.32					± 0.12		
HT HCT	%	37.7	37.7	37.9	38.2	38.3					± 1.0		
	l/l	0.377	0.377	0.379	0.382	0.383					± 0.010		
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	243	241	245	245					± 10		
VPM MPV	µm ³ ; fl	8.8	8.5	N/A	9.2	8.1					± 0.5		

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal


 (Exp.) 2019-03-05
 (YYYY - MM - DD)

LOT		CX 431		CAL		(Exp.) 2019-03-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9		8.9	8.9	9.2	8.8	9.2	8.9	8.9	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50		4.50	4.52	4.54	4.55	4.54	4.51	4.52	± 0.06	
HB HGB	g/dl	13.2		13.3	13.4	13.4	13.3	13.4	13.3	13.4	± 0.2	
	g/l	132		133	134	134	133	134	133	134	± 2	
	mmol/l	8.20		8.26	8.32	8.32	8.26	8.32	8.26	8.32	± 0.12	
HT HCT	%	37.7		36.1	36.7	36.9	36.3	36.7	36.5	36.6	± 1.0	
	l/l	0.377		0.361	0.367	0.369	0.363	0.367	0.365	0.366	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		241	237	241	252	241	250	249	± 10	
VPM MPV	µm ³ ; fl	8.8		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300060880
 BACK / VERSO
 Ref: TEMP-0387 Rev.46

ABX Minocal


LOT
CX 431
CAL

(Exp.) 2019-03-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56									± 0.06
HB HGB	g/dl	13.3									± 0.2
	g/l	133									± 2
	mmol/l	8.26									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	257									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	81.0									± 4.0

1300060884

Ref: TEMP-0387 Rev.46