

ABX Minocal


LOT
CX 459
CAL

(Exp.) 2021-07-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.8			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.56			4.56					± 0.06
HB HGB	g/dl	13.6	13.8			13.7					± 0.2
	g/l	136	138			137					± 2
	mmol/l	8.45	8.57			8.51					± 0.12
HT HCT	%	38.8	39.0			39.0					± 1.0
	l/l	0.388	0.390			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	257	254			254					± 10
VPM MPV	µm ³ ; fl	8.3	8.3			7.8					± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal


 (Exp.) 2021-07-05
 (YYYY - MM - DD)

LOT		CX 459		CAL		<div><div><div></div></div><div>(Exp.) 2021-07-05</div><div>(YYYY - MM - DD)</div></div>					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1		9.0	8.9	9.2	9.0	9.0	9.0	8.9	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55		4.45	4.51	4.55	4.52	4.52	4.55	4.51	± 0.06
HB HGB	g/dl	13.5		13.5	13.6	13.6	13.5	13.6	13.7	13.6	± 0.2
	g/l	135		135	136	136	135	136	137	136	± 2
	mmol/l	8.38		8.38	8.45	8.45	8.38	8.45	8.51	8.45	± 0.12
HT HCT	%	38.4		36.0	37.0	37.1	36.4	36.8	37.3	37.0	± 1.0
	l/l	0.384		0.360	0.370	0.371	0.364	0.368	0.373	0.370	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255		248	243	255	255	255	254	257	± 10
VPM MPV	µm ³ ; fl	8.7		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300060880
 BACK / VERSO
 Ref: TEMP-0387 Rev.46

ABX Minocal



LOT		CX 459		CAL		<div><div></div><div>(Exp.) 2021-07-05</div><div>(YYYY - MM - DD)</div></div>					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN									TOLERANCES TOLERANCE
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58									± 0.06
HB HGB	g/dl	13.6									± 0.2
	g/l	136									± 2
	mmol/l	8.45									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	259									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	82.5									± 4.0

 1300060884
 Ref: TEMP-0387 Rev.46