

ABX Minocal


 (Exp.) 2019-09-05
 (YYYY - MM - DD)

LOT		CX 437		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.5	8.8	8.8	8.9						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.48	4.45	4.45	4.46						± 0.06
HB HGB	g/dl	13.5	13.7	13.5	13.6	13.6						± 0.2
	g/l	135	137	135	136	136						± 2
	mmol/l	8.38	8.51	8.38	8.45	8.45						± 0.12
HT HCT	%	38.3	38.8	38.5	38.7	38.6						± 1.0
	l/l	0.383	0.388	0.385	0.387	0.386						± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	242	235	247	244						± 10
VPM MPV	µm ³ ; fl	8.0	8.0	N/A	8.6	7.6						± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal


 (Exp.) 2019-09-05
 (YYYY - MM - DD)

LOT		CX 437		CAL		(Exp.) 2019-09-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9		9.1	8.9	9.3	8.9	9.3	8.9	8.9	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45		4.40	4.46	4.49	4.44	4.49	4.47	4.46	± 0.06	
HB HGB	g/dl	13.4		13.5	13.7	13.5	13.5	13.5	13.6	13.7	± 0.2	
	g/l	134		135	137	135	135	135	136	137	± 2	
	mmol/l	8.32		8.38	8.51	8.38	8.38	8.38	8.45	8.51	± 0.12	
HT HCT	%	38.3		36.1	37.5	37.3	36.6	37.3	37.3	37.5	± 1.0	
	l/l	0.383		0.361	0.375	0.373	0.366	0.373	0.373	0.375	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242		233	240	248	243	248	240	240	± 10	
VPM MPV	µm ³ ; fl	8.5		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300060880
 BACK / VERSO
 Ref: TEMP-0387 Rev.46

ABX Minocal


LOT
CX 437
CAL

(Exp.) 2019-09-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN									
		H1500 H2500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65									± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51									± 0.06
HB HGB	g/dl	13.5									± 0.2
	g/l	135									± 2
	mmol/l	8.38									± 0.12
HT HCT	%	N/A									N/A
	l/l	N/A									N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245									± 10
VPM MPV	µm ³ ; fl	N/A									N/A
VGM MCV	µm ³ ; fl	83.8									± 4.0

 1300060884
 Ref: TEMP-0387 Rev.46