

ABX Minocal



(Exp.) 2019-08-05

(YYYY - MM - DD)

LOT		CX 436		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		MICROS		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.4	8.7	8.7	8.8						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.55	4.54	4.54	4.57						± 0.06
HB HGB	g/dl	13.5	13.8	13.6	13.7	13.8						± 0.2
	g/l	135	138	136	137	138						± 2
	mmol/l	8.38	8.57	8.45	8.51	8.57						± 0.12
HT HCT	%	38.4	38.7	38.6	38.8	38.8						± 1.0
	l/l	0.384	0.387	0.386	0.388	0.388						± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	237	230	238	236						± 10
VPM MPV	µm ³ ; fl	8.7	8.4	N/A	9.2	8.0						± 0.5

 1300060880
 FRONT / RECTO
 Ref. TEMP-0387 Rev.46

ABX Minocal


 (Exp.) 2019-08-05
 (YYYY - MM - DD)

LOT		CX 436		CAL		(Exp.) 2019-08-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8		8.8	8.8	9.2	8.8	9.2	8.8	8.8	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53		4.53	4.58	4.60	4.55	4.60	4.57	4.58	± 0.06	
HB HGB	g/dl	13.5		13.6	13.8	13.7	13.6	13.7	13.7	13.8	± 0.2	
	g/l	135		136	138	137	136	137	137	138	± 2	
	mmol/l	8.38		8.45	8.57	8.51	8.45	8.51	8.51	8.57	± 0.12	
HT HCT	%	38.1		36.9	37.6	37.7	36.9	37.5	37.3	37.6	± 1.0	
	l/l	0.381		0.369	0.376	0.377	0.369	0.375	0.373	0.376	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240		227	226	240	235	240	235	233	± 10	
VPM MPV	µm ³ ; fl	9.0		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5	

 1300060880
 BACK / VERSO
 Ref. TEMP-0387 Rev.46

ABX Minocal



LOT		CX 436		CAL		(Exp.) 2019-08-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio								TOLERANCES TOLERANCE			
				YUMIZEN											
				H1500 H2500											
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.55										± 0.20		
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63										± 0.06		
HB	HGB	g/dl	13.7										± 0.2		
		g/l	137										± 2		
		mmol/l	8.51										± 0.12		
HT	HCT	%	N/A										N/A		
		l/l	N/A										N/A		
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	246										± 10		
VPM	MPV	µm ³ ; fl	N/A										N/A		
VGM	MCV	µm ³ ; fl	82.0										± 4.0		

 1300060884
 Ref: TEMP-0387 Rev.46