

ABX Minocal



(Exp.) 2021-12-05

(YYYY - MM - DD)

LOT		CX 464		CAL		(Exp.) 2021-12-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.55			4.50					± 0.06
HB HGB	g/dl	13.5	13.6			13.4					± 0.2
	g/l	135	136			134					± 2
	mmol/l	8.38	8.45			8.32					± 0.12
HT HCT	%	38.6	38.7			38.5					± 1.0
	l/l	0.386	0.387			0.385					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	248			245					± 10
VPM MPV	µm ³ ; fl	8.8	8.6			8.1					± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.47

ABX Minocal


 (Exp.) 2021-12-05
 (YYYY - MM - DD)

LOT		CX 464		CAL		(Exp.) 2021-12-05 (YYYY - MM - DD)																	
Rev 1		ABX Lysebio									TOLERANCES TOLERANCE												
PARAMETRES PARAMETERS		MICROS		MICROS		PENTRA		PENTRA		PENTRA													
		EMI CRP		LC-767G		MS CRP		XLR		DX NEXUS DF NEXUS		MS60		DX120 - DF120		60 60 C+ ES60		XL80					
GB WBC		10 ³ /mm ³ ;10 ⁹ /l		9.1		9.0		9.1		9.0		9.4		9.0		9.2		9.0		9.0		± 0.2	
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l		4.46		4.45		4.45		4.48		4.48		4.48		4.51		4.52		4.48		± 0.06	
HB HGB		g/dl		13.3		13.2		13.3		13.4		13.3		13.3		13.4		13.4		13.4		± 0.2	
		g/l		133		132		133		134		133		133		134		134		134		± 2	
		mmol/l		8.26		8.20		8.26		8.32		8.26		8.26		8.32		8.32		8.32		± 0.12	
HT HCT		%		37.7		37.4		36.3		36.5		36.3		36.1		36.5		36.6		36.5		± 1.0	
		l/l		0.377		0.374		0.363		0.365		0.363		0.361		0.365		0.366		0.365		± 0.010	
PLAQ. PLTS		10 ³ /mm ³ ;10 ⁹ /l		249		248		244		240		246		246		250		245		240		± 10	
VPM MPV		µm ³ ; fl		9.1		9.0		N/A		N/A		N/A		N/A		N/A		N/A		N/A		± 0.5	

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300090308

ABX Minocal

LOT		CX 464		CAL		(Exp.) 2021-12-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52									± 0.06	
HB	HGB	g/dl	13.4									± 0.2	
		g/l	134									± 2	
		mmol/l	8.32									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	254									± 10	
VPM	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	82.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.47