

ABX Minocal



LOT

CX 444

CAL



(Exp.) 2020-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP 200	8P	18P	60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.8	8.6	8.8	8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52	4.54	4.50	4.53					± 0.06
HB HGB	g/dl	14.2	14.2	14.1	14.1	14.2					± 0.2
	g/l	142	142	141	141	142					± 2
	mmol/l	8.82	8.82	8.76	8.76	8.82					± 0.12
HT HCT	%	40.2	40.2	40.4	40.1	40.3					± 1.0
	l/l	0.402	0.402	0.404	0.401	0.403					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	242	248	246	246					± 10
VPM MPV	µm ³ ; fl	7.8	7.9	N/A	8.4	7.5					± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060880

ABX Minocal



LOT

CX 444
Rev 1

CAL

(Exp.) **2020-04-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1		9.1	8.9	9.2	9.0	9.2	8.9	8.9	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52		4.52	4.52	4.56	4.50	4.56	4.55	4.52	± 0.06
HB HGB	g/dl	14.0		14.3	14.1	14.2	14.0	14.2	14.2	14.1	± 0.2
	g/l	140		143	141	142	140	142	142	141	± 2
	mmol/l	8.69		8.88	8.76	8.82	8.69	8.82	8.82	8.76	± 0.12
HT HCT	%	39.8		38.4	38.4	39.2	38.0	39.0	38.8	38.6	± 1.0
	l/l	0.398		0.384	0.384	0.392	0.380	0.390	0.388	0.386	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	246		240	247	249	248	249	250	247	± 10
VPM MPV	µm ³ ; fl	8.2		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060880

ABX Minocal



LOT

CX 444

CAL



(Exp.) 2020-04-05
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN								
		H1500 H2500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.65								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56								± 0.06
HB HGB	g/dl	14.2								± 0.2
	g/l	142								± 2
	mmol/l	8.82								± 0.12
HT HCT	%	N/A								N/A
	l/l	N/A								N/A
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248								± 10
VPM MPV	µm ³ ; fl	N/A								N/A
VGM MCV	µm ³ ; fl	85.5								± 4.0

1300060884
Ref: TEMP-0387 Rev.46