

# ABX Minocal



**LOT**

**CX 419**

**CAL**



**(Exp.) 2018-03-05**

Rev 1

(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lyse        |                |              |              |                             |        |        |        | TOLERANCES<br>TOLERANCE |         |
|--------------------------|---|-----------------|----------------|--------------|--------------|-----------------------------|--------|--------|--------|-------------------------|---------|
|                          |   | MICROS          | MICROS         | MICROS       | MICROS       | MICROS                      | PENTRA | PENTRA | PENTRA |                         | PENTRA  |
|                          |   | ES60<br>Care ST | CRP<br>CRP 200 | 8P           | 18P          | 60 8P<br>60 18P<br>ADVIA 60 |        |        |        |                         |         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>9.4</b>      | <b>9.4</b>     | <b>9.1</b>   | <b>9.3</b>   | <b>9.4</b>                  |        |        |        |                         | ± 0.2   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | <b>4.52</b>     | <b>4.54</b>    | <b>4.52</b>  | <b>4.53</b>  | <b>4.52</b>                 |        |        |        |                         | ± 0.06  |
| HB HGB                   | g/dl  | <b>13.8</b>     | <b>13.9</b>    | <b>13.8</b>  | <b>13.8</b>  | <b>13.9</b>                 |        |        |        |                         | ± 0.2   |
|                          | g/l   | <b>138</b>      | <b>139</b>     | <b>138</b>   | <b>138</b>   | <b>139</b>                  |        |        |        |                         | ± 2     |
|                          | mmol/l  | <b>8.57</b>     | <b>8.63</b>    | <b>8.57</b>  | <b>8.57</b>  | <b>8.63</b>                 |        |        |        |                         | ± 0.12  |
| HT HCT                   | %   | <b>39.2</b>     | <b>39.5</b>    | <b>39.3</b>  | <b>39.4</b>  | <b>39.3</b>                 |        |        |        |                         | ± 1.0   |
|                          | l/l   | <b>0.392</b>    | <b>0.395</b>   | <b>0.393</b> | <b>0.394</b> | <b>0.393</b>                |        |        |        |                         | ± 0.010 |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>254</b>      | <b>250</b>     | <b>247</b>   | <b>254</b>   | <b>254</b>                  |        |        |        |                         | ± 10    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | <b>9.0</b>      | <b>8.8</b>     | N/A          | <b>9.2</b>   | <b>8.3</b>                  |        |        |        |                         | ± 0.5   |

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1219930090 -A

# ABX Minocal



**LOT**

**CX 419**

**CAL**



**(Exp.) 2018-03-05**  
(YYYY - MM - DD)

Rev 1

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio  |        |              |              |                      |              |                                  |                     |              | TOLERANCES<br>TOLERANCE |
|--------------------------|---|--------------|--------|--------------|--------------|----------------------|--------------|----------------------------------|---------------------|--------------|-------------------------|
|                          |   | MICROS       | MICROS | PENTRA       | PENTRA       | PENTRA               | PENTRA       | PENTRA                           | PENTRA              | PENTRA       |                         |
|                          |   | EMI CRP      |        | MS CRP       | XLR          | DX NEXUS<br>DF NEXUS | MS60         | 120 - 120 RETIC<br>DX120 - DF120 | 60<br>60 C+<br>ES60 | 80<br>XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>9.5</b>   |        | <b>9.6</b>   | <b>9.2</b>   | <b>9.4</b>           | <b>9.2</b>   | <b>9.4</b>                       | <b>9.3</b>          | <b>9.3</b>   | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | <b>4.51</b>  |        | <b>4.52</b>  | <b>4.53</b>  | <b>4.56</b>          | <b>4.51</b>  | <b>4.56</b>                      | <b>4.52</b>         | <b>4.53</b>  | ± 0.06                  |
| HB HGB                   | g/dl  | <b>13.6</b>  |        | <b>13.8</b>  | <b>13.8</b>  | <b>13.7</b>          | <b>13.8</b>  | <b>13.7</b>                      | <b>13.8</b>         | <b>13.8</b>  | ± 0.2                   |
|                          | g/l   | <b>136</b>   |        | <b>138</b>   | <b>138</b>   | <b>137</b>           | <b>138</b>   | <b>137</b>                       | <b>138</b>          | <b>138</b>   | ± 2                     |
|                          | mmol/l  | <b>8.45</b>  |        | <b>8.57</b>  | <b>8.57</b>  | <b>8.51</b>          | <b>8.57</b>  | <b>8.51</b>                      | <b>8.57</b>         | <b>8.57</b>  | ± 0.12                  |
| HT HCT                   | %   | <b>38.8</b>  |        | <b>37.2</b>  | <b>37.6</b>  | <b>38.3</b>          | <b>37.2</b>  | <b>38.1</b>                      | <b>37.3</b>         | <b>37.7</b>  | ± 1.0                   |
|                          | l/l   | <b>0.388</b> |        | <b>0.372</b> | <b>0.376</b> | <b>0.383</b>         | <b>0.372</b> | <b>0.381</b>                     | <b>0.373</b>        | <b>0.377</b> | ± 0.010                 |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>256</b>   |        | <b>248</b>   | <b>245</b>   | <b>254</b>           | <b>257</b>   | <b>254</b>                       | <b>255</b>          | <b>253</b>   | ± 10                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | <b>9.0</b>   |        | N/A          | N/A          | N/A                  | N/A          | N/A                              | N/A                 | N/A          | ± 0.5                   |

Ref: TEMP-0387 Rev.44 BACK / VERSO 1219930090 -A

# ABX Minocal



**LOT**

**CX 419**

**CAL**



**(Exp.) 2018-03-05**  
(YYYY-MM-DD)

Rev **1**

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio    |  |  |  |  |  |  |  | TOLERANCES<br>TOLERANCE |
|--------------------------|---|----------------|--|--|--|--|--|--|--|-------------------------|
|                          |   | YUMIZEN        |  |  |  |  |  |  |  |                         |
|                          |   | H1500<br>H2500 |  |  |  |  |  |  |  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>9.1</b>     |  |  |  |  |  |  |  | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | <b>4.60</b>    |  |  |  |  |  |  |  | ± 0.06                  |
| HB HGB                   | g/dl  | <b>13.8</b>    |  |  |  |  |  |  |  | ± 0.2                   |
|                          | g/l   | <b>138</b>     |  |  |  |  |  |  |  | ± 2                     |
|                          | mmol/l  | <b>8.57</b>    |  |  |  |  |  |  |  | ± 0.12                  |
| HT HCT                   | %   | N/A            |  |  |  |  |  |  |  | N/A                     |
|                          | l/l   | N/A            |  |  |  |  |  |  |  | N/A                     |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>267</b>     |  |  |  |  |  |  |  | ± 10                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | N/A            |  |  |  |  |  |  |  | N/A                     |
| VGM MCV                  | µm <sup>3</sup> ; fl                                  | <b>83.0</b>    |  |  |  |  |  |  |  | ± 4                     |

1300032517-A  
Ref: TEMP-0387 Rev.44