

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT CX 053
Rev 1

CAL

 (Exp.) 2013-06-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse										TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.5	9.0	9.1	9.0	9.2		9.5	9.3	9.3	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.51	4.53	4.53	4.51		4.54	4.53	4.54	± 0.06	
HB HGB	g/dl	13.6	13.5	13.3	13.5	13.5		13.4	13.4	13.5	± 0.2	
	g/l	136	135	133	135	135		134	134	135	± 2	
	mmol/l	8.45	8.38	8.26	8.38	8.38		8.32	8.32	8.38	± 0.12	
HT HCT	%	38.7	38.3	38.7	38.7	38.6		37.2	36.9	37.3	± 1.0	
	l/l	0.387	0.383	0.387	0.387	0.386		0.372	0.369	0.373	± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	251	240	254	253	245		248	247	248	± 10	
VPM MPV	µm ³ ; fl	7.8	7.8	N/A	7.8	7.3		N/A	N/A	N/A	± 0.5	

Ref: TEMP-0387 Rev.32 FRONT / RECTO 9930081 -A

ABX Minocal



LOT

CX 053

CAL



(Exp.) 2013-06-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	CRP CRP 200	8P	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4			9.3	9.5	9.3	9.5	9.3	9.3	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51			4.54	4.54	4.50	4.54	4.53	4.54	± 0.06
HB HGB	g/dl	13.4			13.5	13.4	13.4	13.4	13.4	13.5	± 0.2
	g/l	134			135	134	134	134	134	135	± 2
	mmol/l	8.32			8.38	8.32	8.32	8.32	8.32	8.38	± 0.12
HT HCT	%	38.1			37.3	37.2	36.0	37.2	36.9	37.3	± 1.0
	l/l	0.381			0.373	0.372	0.360	0.372	0.369	0.373	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248			248	248	250	248	247	248	± 10
VPM MPV	µm ³ ; fl	7.9			N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.32 BACK / VERSO 9930081 -A