

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT CX 043
Rev 1

CAL

 (Exp.) 2013-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.7	9.5	9.5	9.4	9.5		9.7	9.5	9.5	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.53	4.51	4.51	4.53		4.56	4.53	4.54	± 0.06
HB HGB	g/dl	13.7	13.6	13.6	13.6	13.6		13.6	13.6	13.6	± 0.2
	g/l	137	136	136	136	136		136	136	136	± 2
	mmol/l	8.51	8.45	8.45	8.45	8.45		8.45	8.45	8.45	± 0.12
HT HCT	%	38.8	39.0	39.2	39.0	39.0		37.8	37.4	37.6	± 1.0
	l/l	0.388	0.390	0.392	0.390	0.390		0.378	0.374	0.376	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248	250	250	254	245		250	250	250	± 10
VPM MPV	µm ³ ; fl	7.8	7.9	N/A	8.1	7.5		N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.31 FRONT / RECTO 9930077 -A

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT CX 043
Rev 1

CAL

 (Exp.) 2013-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	CRP CRP 200	8P	18P	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3				9.7	9.4	9.7	9.5	9.5	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58				4.56	4.51	4.56	4.53	4.54	± 0.06
HB HGB	g/dl	13.4				13.6	13.6	13.6	13.6	13.6	± 0.2
	g/l	134				136	136	136	136	136	± 2
	mmol/l	8.32				8.45	8.45	8.45	8.45	8.45	± 0.12
HT HCT	%	39.2				37.8	36.4	37.8	37.4	37.6	± 1.0
	l/l	0.392				0.378	0.364	0.378	0.374	0.376	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238				250	250	250	250	250	± 10
VPM MPV	µm ³ ; fl	8.0				N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.31 BACK / VERSO 9930077 -A