

# ABX Minocal

| LOT                      |                                                       | CX 417          |                | CAL    |        | <div><div></div><div>(Exp.) 2018-01-05</div><div>(YYYY - MM - DD)</div></div> |        |        |        |        |                         |
|--------------------------|-------------------------------------------------------|-----------------|----------------|--------|--------|-------------------------------------------------------------------------------|--------|--------|--------|--------|-------------------------|
| Rev 1                    |                                                       |                 |                |        |        |                                                                               |        |        |        |        |                         |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lyse        |                |        |        |                                                                               |        |        |        |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS          | MICROS         | MICROS | MICROS | MICROS                                                                        | PENTRA | PENTRA | PENTRA | PENTRA |                         |
|                          |                                                       | ES60<br>Care ST | CRP<br>CRP 200 | 8P     | 18P    | 60 8P<br>60 18P<br>ADVIA 60                                                   |        |        |        |        |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.3             | 9.3            | 9.0    | 9.3    | 9.3                                                                           |        |        |        |        | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.67            | 4.63           | 4.70   | 4.63   | 4.65                                                                          |        |        |        |        | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.7            | 13.7           | 13.7   | 13.6   | 13.7                                                                          |        |        |        |        | ± 0.2                   |
|                          | g/l                                                   | 137             | 137            | 137    | 136    | 137                                                                           |        |        |        |        | ± 2                     |
|                          | mmol/l                                                | 8.51            | 8.51           | 8.51   | 8.45   | 8.51                                                                          |        |        |        |        | ± 0.12                  |
| HT HCT                   | %                                                     | 39.0            | 38.7           | 39.2   | 38.7   | 39.0                                                                          |        |        |        |        | ± 1.0                   |
|                          | l/l                                                   | 0.390           | 0.387          | 0.392  | 0.387  | 0.390                                                                         |        |        |        |        | ± 0.010                 |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 254             | 251            | 253    | 253    | 253                                                                           |        |        |        |        | ± 10                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | 8.8             | 8.6            | N/A    | 8.8    | 8.1                                                                           |        |        |        |        | ± 0.5                   |

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1219930090 -A

# ABX Minocal


 (Exp.) 2018-01-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 417      |        | CAL    |        | <div><div><div></div></div><div>(Exp.) 2018-01-05</div><div>(YYYY - MM - DD )</div></div> |        |                                  |                     |            |                         |  |  |
|--------------------------|-------------------------------------------------------|-------------|--------|--------|--------|-------------------------------------------------------------------------------------------|--------|----------------------------------|---------------------|------------|-------------------------|--|--|
|                          |                                                       | Rev 1       |        |        |        |                                                                                           |        |                                  |                     |            |                         |  |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |        |        |        |                                                                                           |        |                                  |                     |            | TOLERANCES<br>TOLERANCE |  |  |
|                          |                                                       | MICROS      | MICROS | PENTRA | PENTRA | PENTRA                                                                                    | PENTRA | PENTRA                           | PENTRA              | PENTRA     |                         |  |  |
|                          |                                                       | EMI CRP     |        | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                                      | MS60   | 120 - 120 RETIC<br>DX120 - DF120 | 60<br>60 C+<br>ES60 | 80<br>XL80 |                         |  |  |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.3         |        | 9.6    | 9.3    | 9.5                                                                                       | 9.3    | 9.5                              | 9.3                 | 9.3        | ± 0.2                   |  |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.65        |        | 4.66   | 4.66   | 4.70                                                                                      | 4.72   | 4.70                             | 4.67                | 4.68       | ± 0.06                  |  |  |
| HB HGB                   | g/dl                                                  | 13.5        |        | 13.6   | 13.6   | 13.7                                                                                      | 13.8   | 13.7                             | 13.6                | 13.6       | ± 0.2                   |  |  |
|                          | g/l                                                   | 135         |        | 136    | 136    | 137                                                                                       | 138    | 137                              | 136                 | 136        | ± 2                     |  |  |
|                          | mmol/l                                                | 8.38        |        | 8.45   | 8.45   | 8.51                                                                                      | 8.57   | 8.51                             | 8.45                | 8.45       | ± 0.12                  |  |  |
| HT HCT                   | %                                                     | 38.5        |        | 36.8   | 37.0   | 38.1                                                                                      | 37.3   | 37.8                             | 37.1                | 37.4       | ± 1.0                   |  |  |
|                          | l/l                                                   | 0.385       |        | 0.368  | 0.370  | 0.381                                                                                     | 0.373  | 0.378                            | 0.371               | 0.374      | ± 0.010                 |  |  |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 257         |        | 250    | 250    | 255                                                                                       | 261    | 255                              | 258                 | 255        | ± 10                    |  |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | 8.7         |        | N/A    | N/A    | N/A                                                                                       | N/A    | N/A                              | N/A                 | N/A        | ± 0.5                   |  |  |

 1219930090 -A  
 BACK / VERSO  
 Ref: TEMP-0387 Rev.44

# ABX Minocal



| LOT                      |                                                       | CX 417         |  | CAL |  | <div><div></div><div>(Exp.) 2018-01-05</div><div>(YYYY - MM - DD)</div></div> |  |  |  |  |                         |  |  |
|--------------------------|-------------------------------------------------------|----------------|--|-----|--|-------------------------------------------------------------------------------|--|--|--|--|-------------------------|--|--|
|                          |                                                       | Rev 1          |  |     |  |                                                                               |  |  |  |  |                         |  |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio    |  |     |  |                                                                               |  |  |  |  | TOLERANCES<br>TOLERANCE |  |  |
|                          |                                                       | YUMIZEN        |  |     |  |                                                                               |  |  |  |  |                         |  |  |
|                          |                                                       | H1500<br>H2500 |  |     |  |                                                                               |  |  |  |  |                         |  |  |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.1            |  |     |  |                                                                               |  |  |  |  | ± 0.2                   |  |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.74           |  |     |  |                                                                               |  |  |  |  | ± 0.06                  |  |  |
| HB HGB                   | g/dl                                                  | 13.6           |  |     |  |                                                                               |  |  |  |  | ± 0.2                   |  |  |
|                          | g/l                                                   | 136            |  |     |  |                                                                               |  |  |  |  | ± 2                     |  |  |
|                          | mmol/l                                                | 8.45           |  |     |  |                                                                               |  |  |  |  | ± 0.12                  |  |  |
| HT HCT                   | %                                                     | N/A            |  |     |  |                                                                               |  |  |  |  | N/A                     |  |  |
|                          | l/l                                                   | N/A            |  |     |  |                                                                               |  |  |  |  | N/A                     |  |  |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 270            |  |     |  |                                                                               |  |  |  |  | ± 10                    |  |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | N/A            |  |     |  |                                                                               |  |  |  |  | N/A                     |  |  |
| VGM MCV                  | µm <sup>3</sup> ; fl                                  | 79             |  |     |  |                                                                               |  |  |  |  | ± 4                     |  |  |

 1300032517-A  
 Ref: TEMP-0387 Rev.44