

# ABX Minocal



**LOT** CX 024  
Rev 2

**CAL**

 (Exp.) 2014-03-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60		80 XL80
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	9.5	9.2	9.2	9.2	9.4			9.4	9.4	± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.46	4.40	4.42	4.42	4.44			4.46	4.45	± 0.06
HB HGB	g/dl	13.7	13.6	13.4	13.6	13.6			13.6	13.6	± 0.2
	g/l	137	136	134	136	136			136	136	± 2
	mmol/l	8.51	8.45	8.32	8.45	8.45			8.45	8.45	± 0.12
HT HCT	%	38.7	38.5	38.5	38.8	38.5			37.3	37.4	± 1.0
	l/l	0.387	0.385	0.385	0.388	0.385			0.373	0.374	± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	257	250	248	259	256			258	258	± 10
VPM MPV	µm <sup>3</sup> ; fl	8.0	7.9	N/A	8.0	7.5			N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.35 FRONT / RECTO 9930081 -B

**LOT** CX 024  
Rev 2

**CAL**

 (Exp.) 2014-03-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	CRP CRP 200	8P	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	9.2			9.4	9.5	9.3	9.5	9.4	9.4	± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.43			4.45	4.44	4.42	4.44	4.46	4.45	± 0.06
HB HGB	g/dl	13.5			13.6	13.6	13.6	13.6	13.6	13.6	± 0.2
	g/l	135			136	136	136	136	136	136	± 2
	mmol/l	8.38			8.45	8.45	8.45	8.45	8.45	8.45	± 0.12
HT HCT	%	38.3			37.4	37.3	36.0	37.3	37.3	37.4	± 1.0
	l/l	0.383			0.374	0.373	0.360	0.373	0.373	0.374	± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	250			258	252	257	252	258	258	± 10
VPM MPV	µm <sup>3</sup> ; fl	8.2			N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.35 BACK / VERSO 9930081 -B

**LOT** CX 024  
Rev 2

**CAL**

 (Exp.) 2014-03-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	9.2	9.2			9.2					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.43	4.40			4.44					± 0.06
HB HGB	g/dl	13.5	13.3			13.2					± 0.2
	g/l	135	133			132					± 2
	mmol/l	8.38	8.28			8.19					± 0.12
HT HCT	%	38.3	38.5			38.5					± 1.0
	l/l	0.383	0.385			0.385					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	250	250			256					± 10
VPM MPV	µm <sup>3</sup> ; fl	8.2	7.9			7.5					± 0.5

Ref: TEMP-0387 Rev.35 9930088-B