

ABX Minocal

LOT

CX 405

CAL



(Exp.) 2017-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.5	9.5			9.6					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.61	4.58			4.59					± 0.06
HB HGB	g/dl	13.8	13.5			13.9					± 0.2
	g/l	138	135			139					± 2
	mmol/l	8.57	8.40			8.63					± 0.12
HT HCT	%	39.4	39.4			39.2					± 1.0
	l/l	0.394	0.394			0.392					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	245			248					± 10
VPM MPV	µm ³ ; fl	8.7	8.8			9.1					± 0.5

ABX Minocal

LOT

CX 405

Rev 1

CAL



(Exp.) **2017-01-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.2	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.61	4.60							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.2	39.1							± 1.0
	l/l	0.392	0.391							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	248							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.43 FRONT / RECTO 1300027265-A