

ABX Minocal

LOT

CX 454

CAL



(Exp.) 2021-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.8			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.52			4.55					± 0.06
HB HGB	g/dl	13.9	13.4			13.8					± 0.2
	g/l	139	134			138					± 2
	mmol/l	8.63	8.34			8.57					± 0.12
HT HCT	%	39.0	38.9			39.6					± 1.0
	l/l	0.390	0.389			0.396					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	242			244					± 10
VPM MPV	µm ³ ; fl	8.7	8.2			9.2					± 0.5

ABX Minocal

LOT

CX 454

CAL



(Exp.) 2021-02-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.52							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	39.0	38.9							± 1.0
	l/l	0.390	0.389							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	249							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881