

ABX Minocal

LOT

CX 444

CAL



(Exp.) 2020-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.8			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52			4.55					± 0.06
HB HGB	g/dl	14.3	13.9			14.3					± 0.2
	g/l	143	139			143					± 2
	mmol/l	8.88	8.64			8.88					± 0.12
HT HCT	%	40.2	40.2			40.5					± 1.0
	l/l	0.402	0.402			0.405					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	242			248					± 10
VPM MPV	µm ³ ; fl	8.3	7.9			8.7					± 0.5

ABX Minocal

LOT

CX 444

CAL



(Exp.) 2020-04-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.52							± 0.06
HB HGB	g/dl	14.1	14.0							± 0.2
	g/l	141	140							± 2
	mmol/l	8.76	8.69							± 0.12
HT HCT	%	40.3	40.0							± 1.0
	l/l	0.403	0.400							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	254							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881