

ABX Minocal

LOT

CX 414

CAL



(Exp.) 2017-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.40			4.40					± 0.06
HB HGB	g/dl	13.8	13.3			13.8					± 0.2
	g/l	138	133			138					± 2
	mmol/l	8.57	8.28			8.57					± 0.12
HT HCT	%	38.9	38.7			38.6					± 1.0
	l/l	0.389	0.387			0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248	250			255					± 10
VPM MPV	µm ³ ; fl	8.7	8.5			8.6					± 0.5

ABX Minocal

LOT

CX 414

Rev 1

CAL



(Exp.) 2017-10-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.44							± 0.06
HB HGB	g/dl	13.4	13.5							± 0.2
	g/l	134	135							± 2
	mmol/l	8.32	8.38							± 0.12
HT HCT	%	38.6	38.5							± 1.0
	l/l	0.386	0.385							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	253							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A