

ABX Minocal

LOT
CX 425

 Rev **1**
CAL

(Exp.) 2018-09-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.65	4.67			4.65					± 0.06
HB HGB	g/dl	13.8	13.5			13.8					± 0.2
	g/l	138	135			138					± 2
	mmol/l	8.57	8.40			8.57					± 0.12
HT HCT	%	39.1	39.4			39.1					± 1.0
	l/l	0.391	0.394			0.391					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	245			246					± 10
VPM MPV	µm ³ ; fl	8.6	8.0			8.2					± 0.5

ABX Minocal

LOT
CX 425

 Rev **1**
CAL

(Exp.) 2018-09-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62								± 0.06
HB HGB	g/dl	13.5	13.5								± 0.2
	g/l	135	135								± 2
	mmol/l	8.38	8.38								± 0.12
HT HCT	%	38.8	38.8								± 1.0
	l/l	0.388	0.388								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	226	242								± 10
VPM MPV	µm ³ ; fl	N/A	8.3								± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A