

ABX Minocal

LOT

CX 425

CAL



(Exp.) 2018-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.65	4.67			4.65					± 0.06
HB HGB	g/dl	13.8	13.5			13.8					± 0.2
	g/l	138	135			138					± 2
	mmol/l	8.57	8.40			8.57					± 0.12
HT HCT	%	39.1	39.4			39.1					± 1.0
	l/l	0.391	0.394			0.391					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	245			246					± 10
VPM MPV	µm ³ ; fl	8.6	8.0			8.2					± 0.5

ABX Minocal

LOT

CX 425

Rev 1

CAL



(Exp.) 2018-09-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
HT HCT	%	38.8	38.8							± 1.0
	l/l	0.388	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	226	242							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A