

# ABX Minocal

**LOT**

**CX 433**

**CAL**



**(Exp.) 2019-05-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>8.8</b>			<b>9.1</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.42</b>	<b>4.38</b>			<b>4.40</b>					± 0.06
HB HGB	g/dl	<b>13.4</b>	<b>13.2</b>			<b>13.4</b>					± 0.2
	g/l	<b>134</b>	<b>132</b>			<b>134</b>					± 2
	mmol/l	<b>8.32</b>	<b>8.22</b>			<b>8.32</b>					± 0.12
HT HCT	%	<b>37.8</b>	<b>37.7</b>			<b>37.7</b>					± 1.0
	l/l	<b>0.378</b>	<b>0.377</b>			<b>0.377</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>234</b>	<b>248</b>			<b>249</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>8.7</b>	<b>8.0</b>			<b>8.5</b>					± 0.5

# ABX Minocal

**LOT**

**CX 433**

Rev 1

**CAL**



(Exp.) **2019-05-05**  
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.6</b>	<b>8.7</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.40</b>	<b>4.38</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.1</b>							± 0.2
	g/l	<b>132</b>	<b>131</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.14</b>							± 0.12
HT HCT	%	<b>37.8</b>	<b>37.7</b>							± 1.0
	l/l	<b>0.378</b>	<b>0.377</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>232</b>	<b>242</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.6</b>							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881