

ABX Minocal

LOT

CX 403

CAL



(Exp.) 2016-11-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.0			8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.52			4.50					± 0.06
HB HGB	g/dl		13.4			13.3					± 0.2
	g/l		134			133					± 2
	mmol/l		8.34			8.25					± 0.12
HT HCT	%		38.9			38.7					± 1.0
	l/l		0.389			0.387					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		245			248					± 10
VPM MPV	µm ³ ; fl		8.1			7.3					± 0.5

ABX Minocal

LOT

CX 403

Rev 1

CAL



(Exp.) 2016-11-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.53							± 0.06
HB HGB	g/dl	13.5	13.6							± 0.2
	g/l	135	136							± 2
	mmol/l	8.38	8.45							± 0.12
HT HCT	%	38.8	38.7							± 1.0
	l/l	0.388	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233	250							± 10
VPM MPV	µm ³ ; fl	N/A	7.9							± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A