

ABX Minocal

LOT		CX 452		CAL		(Exp.) 2020-12-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.55			4.55					± 0.06
HB HGB	g/dl	14.0	13.8			13.9					± 0.2
	g/l	140	138			139					± 2
	mmol/l	8.69	8.58			8.63					± 0.12
HT HCT	%	39.4	40.0			40.0					± 1.0
	l/l	0.394	0.400			0.400					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	237			240					± 10
VPM MPV	µm ³ ; fl	8.3	7.9			8.6					± 0.5

ABX Minocal



LOT		CX 452		CAL		(Exp.) 2020-12-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.6	8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.50								± 0.06
HB HGB	g/dl	13.7	13.7								± 0.2
	g/l	137	137								± 2
	mmol/l	8.51	8.51								± 0.12
HT HCT	%	39.5	39.6								± 1.0
	l/l	0.395	0.396								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	238								± 10
VPM MPV	µm ³ ; fl	N/A	8.7								± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881