

ABX Minocal

LOT

CX 422

CAL



(Exp.) 2018-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40			4.38					± 0.06
HB HGB	g/dl	12.9	12.5			12.8					± 0.2
	g/l	129	125			128					± 2
	mmol/l	8.01	7.79			7.95					± 0.12
HT HCT	%	36.5	36.7			36.4					± 1.0
	l/l	0.365	0.367			0.364					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	254			253					± 10
VPM MPV	µm ³ ; fl	8.8	8.4			8.5					± 0.5

ABX Minocal

LOT

CX 422

CAL



(Exp.) 2018-06-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.42							± 0.06
HB HGB	g/dl	12.7	12.6							± 0.2
	g/l	127	126							± 2
	mmol/l	7.89	7.82							± 0.12
HT HCT	%	36.7	36.6							± 1.0
	l/l	0.367	0.366							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	236	252							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A