

ABX Minocal



(Exp.) 2017-11-05
(YYYY - MM - DD)

LOT

CX 415

CAL

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.44			4.44					± 0.06
HB HGB	g/dl	13.9	13.4			13.8					± 0.2
	g/l	139	134			138					± 2
	mmol/l	8.63	8.34			8.57					± 0.12
HT HCT	%	39.3	38.9			39.0					± 1.0
	l/l	0.393	0.389			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	264	255			264					± 10
VPM MPV	µm ³ ; fl	7.9	7.7			7.9					± 0.5

ABX Minocal

LOT

CX 415

Rev 1

CAL



(Exp.) 2017-11-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.48							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.1	39.0							± 1.0
	l/l	0.391	0.390							± 0.010
PLAQ. PLTS	10 ⁹ /mm ³ ;10 ⁹ /l	244	255							± 10
VPM MPV	µm ³ ; fl	N/A	7.8							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A