

ABX Minocal


 (Exp.) 2019-07-05
 (YYYY - MM - DD)

LOT		CX 435		CAL		(Exp.) 2019-07-05 (YYYY - MM - DD)						
		Rev 1		Minilysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.4			8.9					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.61	4.55			4.58					± 0.06	
HB HGB	g/dl	13.6	13.2			13.6					± 0.2	
	g/l	136	132			136					± 2	
	mmol/l	8.45	8.22			8.45					± 0.12	
HT HCT	%	38.6	38.2			38.5					± 1.0	
	l/l	0.386	0.382			0.385					± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	242			245					± 10	
VPM MPV	µm ³ ; fl	8.3	7.9			8.4					± 0.5	

ABX Minocal

LOT		CX 435		CAL		(Exp.) 2019-07-05 (YYYY - MM - DD)					
Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS								TOLERANCES TOLERANCE
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.6	8.6								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.57								± 0.06
HB HGB	g/dl	13.5	13.2								± 0.2
	g/l	135	132								± 2
	mmol/l	8.38	8.20								± 0.12
HT HCT	%	38.6	38.4								± 1.0
	l/l	0.386	0.384								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	245								± 10
VPM MPV	µm ³ ; fl	N/A	8.4								± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881