

ABX Minocal

LOT

CX 470

CAL



(Exp.) 2022-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.50			4.46					± 0.06
HB HGB	g/dl	13.8	13.7			14.0					± 0.2
	g/l	138	137			140					± 2
	mmol/l	8.57	8.52			8.69					± 0.12
HT HCT	%	38.9	39.8			40.6					± 1.0
	l/l	0.389	0.398			0.406					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	232	241			240					± 10
VPM MPV	µm ³ ; fl	9.0	8.4			9.0					± 0.5

ABX Minocal

LOT

CX 470

CAL



(Exp.) 2022-06-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.46							± 0.06
HB HGB	g/dl	13.5	13.6							± 0.2
	g/l	135	136							± 2
	mmol/l	8.38	8.45							± 0.12
HT HCT	%	38.9	39.2							± 1.0
	l/l	0.389	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237	245							± 10
VPM MPV	µm ³ ; fl	N/A	9.2							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881