

ABX Minocal

LOT

CX 469

CAL



(Exp.) 2022-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.9			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.61			4.56					± 0.06
HB HGB	g/dl	13.8	13.6			13.8					± 0.2
	g/l	138	136			138					± 2
	mmol/l	8.57	8.46			8.57					± 0.12
HT HCT	%	39.2	39.6			40.1					± 1.0
	l/l	0.392	0.396			0.401					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	241	248			245					± 10
VPM MPV	µm ³ ; fl	8.6	8.0			8.8					± 0.5

ABX Minocal

LOT

CX 469

CAL



(Exp.) 2022-05-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.60							± 0.06
HB HGB	g/dl	13.4	13.5							± 0.2
	g/l	134	135							± 2
	mmol/l	8.32	8.38							± 0.12
HT HCT	%	38.7	39.1							± 1.0
	l/l	0.387	0.391							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	239	250							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881