

ABX Minocal

LOT

CX 458

CAL



(Exp.) 2021-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.54			4.58					± 0.06
HB HGB	g/dl	14.1	13.8			14.1					± 0.2
	g/l	141	138			141					± 2
	mmol/l	8.76	8.58			8.76					± 0.12
HT HCT	%	39.6	39.5			40.8					± 1.0
	l/l	0.396	0.395			0.408					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	238			241					± 10
VPM MPV	µm ³ ; fl	8.7	8.6			9.8					± 0.5

ABX Minocal

LOT

CX 458

Rev 1

CAL



(Exp.) **2021-06-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.56							± 0.06
HB HGB	g/dl	13.9	13.8							± 0.2
	g/l	139	138							± 2
	mmol/l	8.63	8.57							± 0.12
HT HCT	%	40.0	39.4							± 1.0
	l/l	0.400	0.394							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	243							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881