

ABX Minocal

LOT

CX 447

CAL



(Exp.) 2020-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.2			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.46			4.45					± 0.06
HB HGB	g/dl	13.5	13.1			13.4					± 0.2
	g/l	135	131			134					± 2
	mmol/l	8.38	8.15			8.32					± 0.12
HT HCT	%	38.3	37.9			38.5					± 1.0
	l/l	0.383	0.379			0.385					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	247			245					± 10
VPM MPV	µm ³ ; fl	8.2	7.7			8.8					± 0.5

ABX Minocal

LOT

CX 447

CAL



(Exp.) 2020-07-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.0	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.45							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	37.9	37.8							± 1.0
	l/l	0.379	0.378							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	246							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881