

ABX Minocal

LOT

CX 427

CAL



(Exp.) 2018-11-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.54			4.54					± 0.06
HB HGB	g/dl	13.8	13.4			13.7					± 0.2
	g/l	138	134			137					± 2
	mmol/l	8.57	8.34			8.51					± 0.12
HT HCT	%	38.8	38.6			38.6					± 1.0
	l/l	0.388	0.386			0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237	242			245					± 10
VPM MPV	µm ³ ; fl	9.5	9.1			9.4					± 0.5

ABX Minocal

LOT

CX 427
Rev 1

CAL



(Exp.) 2018-11-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.0	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.56							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	38.8	39.0							± 1.0
	l/l	0.388	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	228	235							± 10
VPM MPV	µm ³ ; fl	N/A	9.6							± 0.5

Ref: TEMP-0387 Rev.45 FRONT / RECTO 1300032436-A