

ABX Minocal

LOT

CX 456

CAL



(Exp.) 2021-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.41			4.48					± 0.06
HB HGB	g/dl	13.7	13.3			13.9					± 0.2
	g/l	137	133			139					± 2
	mmol/l	8.51	8.28			8.63					± 0.12
HT HCT	%	38.5	38.6			39.9					± 1.0
	l/l	0.385	0.386			0.399					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	248			252					± 10
VPM MPV	µm ³ ; fl	8.9	8.2			9.5					± 0.5

ABX Minocal

LOT

CX 456

Rev 1

CAL



(Exp.) **2021-04-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.41							± 0.06
HB HGB	g/dl	13.5	13.3							± 0.2
	g/l	135	133							± 2
	mmol/l	8.38	8.26							± 0.12
HT HCT	%	39.2	38.4							± 1.0
	l/l	0.392	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	255							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881