

ABX Minocal

LOT

CX 446

CAL



(Exp.) 2020-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.50			4.54					± 0.06
HB HGB	g/dl	13.2	12.8			13.3					± 0.2
	g/l	132	128			133					± 2
	mmol/l	8.20	7.97			8.26					± 0.12
HT HCT	%	37.4	37.1			37.7					± 1.0
	l/l	0.374	0.371			0.377					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	244			242					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			9.7					± 0.5

ABX Minocal

LOT

CX 446

CAL



(Exp.) 2020-06-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.53							± 0.06
HB HGB	g/dl	13.0	12.9							± 0.2
	g/l	130	129							± 2
	mmol/l	8.07	8.01							± 0.12
HT HCT	%	37.5	37.4							± 1.0
	l/l	0.375	0.374							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	243							± 10
VPM MPV	µm ³ ; fl	N/A	9.3							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881