

ABX Minocal

LOT

CX 416

CAL



(Exp.) 2017-12-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.79	4.79			4.77					± 0.06
HB HGB	g/dl	14.0	13.6			14.1					± 0.2
	g/l	140	136			141					± 2
	mmol/l	8.69	8.46			8.76					± 0.12
HT HCT	%	39.9	40.0			39.8					± 1.0
	l/l	0.399	0.400			0.398					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	251			260					± 10
VPM MPV	µm ³ ; fl	9.5	9.2			9.4					± 0.5

ABX Minocal

LOT

CX 416

Rev 1

CAL



(Exp.) 2017-12-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.79	4.84							± 0.06
HB HGB	g/dl	13.7	13.8							± 0.2
	g/l	137	138							± 2
	mmol/l	8.51	8.57							± 0.12
HT HCT	%	39.9	39.9							± 1.0
	l/l	0.399	0.399							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	250							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A