

ABX Minocal

LOT

CX 105

CAL



(Exp.) 2015-11-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.46			4.46					± 0.06
HB HGB	g/dl		13.6			13.5					± 0.2
	g/l		136			135					± 2
	mmol/l		8.46			8.37					± 0.12
HT HCT	%		38.8			38.8					± 1.0
	l/l		0.388			0.388					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		248			250					± 10
VPM MPV	µm ³ ; fl		8.3			7.7					± 0.5

ABX Minocal

LOT

CX 105

Rev 1

CAL



(Exp.) **2015-11-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.0	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.49							± 0.06
HB HGB	g/dl	13.9	13.6							± 0.2
	g/l	139	136							± 2
	mmol/l	8.63	8.45							± 0.12
HT HCT	%	39.3	39.0							± 1.0
	l/l	0.393	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	247	253							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A