

ABX Minocal

LOT

CX 085

CAL



(Exp.) 2015-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.4			9.3					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.53			4.56					± 0.06
HB HGB	g/dl		13.4			13.4					± 0.2
	g/l		134			134					± 2
	mmol/l		8.34			8.31					± 0.12
HT HCT	%		38.9			39.2					± 1.0
	l/l		0.389			0.392					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		245			248					± 10
VPM MPV	µm ³ ; fl		8.1			7.7					± 0.5

ABX Minocal

LOT

CX 085

Rev 1

CAL



(Exp.) **2015-09-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.2	9.3							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.58							± 0.06
HB HGB	g/dl	13.8	13.6							± 0.2
	g/l	138	136							± 2
	mmol/l	8.57	8.45							± 0.12
HT HCT	%	39.4	38.9							± 1.0
	l/l	0.394	0.389							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	238	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A