

ABX Minocal

LOT		CX 075		CAL		 (Exp.) 2015-08-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.48			4.49					± 0.06
HB HGB	g/dl		13.6			13.5					± 0.2
	g/l		136			135					± 2
	mmol/l		8.46			8.37					± 0.12
HT HCT	%		39.0			39.0					± 1.0
	l/l		0.390			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		242			240					± 10
VPM MPV	µm ³ ; fl		8.4			7.9					± 0.5

ABX Minocal



LOT		CX 075		CAL		(Exp.) 2015-08-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.52								± 0.06
HB HGB	g/dl	13.8	13.7								± 0.2
	g/l	138	137								± 2
	mmol/l	8.57	8.51								± 0.12
HT HCT	%	39.4	39.1								± 1.0
	l/l	0.394	0.391								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	232	245								± 10
VPM MPV	µm ³ ; fl	N/A	8.5								± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A