

ABX Minocal

LOT

CX 441

CAL



(Exp.) 2020-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.3			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.48			4.51					± 0.06
HB HGB	g/dl	14.0	13.7			14.1					± 0.2
	g/l	140	137			141					± 2
	mmol/l	8.69	8.52			8.76					± 0.12
HT HCT	%	39.5	39.2			39.6					± 1.0
	l/l	0.395	0.392			0.396					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	240			240					± 10
VPM MPV	µm ³ ; fl	9.3	8.9			9.4					± 0.5

ABX Minocal

LOT

CX 441
Rev 1

CAL



(Exp.) **2020-01-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.52							± 0.06
HB HGB	g/dl	13.9	13.8							± 0.2
	g/l	139	138							± 2
	mmol/l	8.63	8.57							± 0.12
HT HCT	%	39.5	39.6							± 1.0
	l/l	0.395	0.396							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	225	243							± 10
VPM MPV	µm ³ ; fl	N/A	9.7							± 0.5

Ref: TEMP-0387 Rev.46 FRONT / RECTO 1300060881